

Antibodies Product Datasheet

Product Name:	PSMC3 Antibody
Catalog Number:	C14013
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PRS6A. AA range:271-320
Concentration:	1 mg/ml
Gene Name:	PSMC3
Protein Name:	26S protease regulatory subunit 6A
Other Name:	PSMC3; TBP1; 26S protease regulatory subunit 6A; 26S proteasome AAA-ATPase subunit RPT5; Proteasome 26S subunit ATPase 3; Proteasome subunit P50; Tat-binding protein 1; TBP-1
Human Swiss Prot No.:	P17980
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PSMC3 Polyclonal Antibody detects endogenous levels of PSMC3 protein.
Storage and Stability:	-20°C/1 year