

Antibodies Product Datasheet

Product Name:	RUNX3 Antibody
Catalog Number:	R12-2341
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RUNX3. AA range:133-182
Concentration:	1 mg/ml
Gene Name:	RUNX3
Protein Name:	Runt-related transcription factor 3
Other Name:	RUNX3; AML2; CBFA3; PEBP2A3; Runt-related transcription factor 3; Acute myeloid leukemia 2 protein; Core-binding factor subunit alpha-3; CBF-alpha-3; Oncogene AML-2; Polyomavirus enhancer-binding protein 2 alpha C subunit; PEA2-alpha C; PEB
Human Swiss Prot No.:	Q13761
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RUNX3 Polyclonal Antibody detects endogenous levels of RUNX3 protein.
Storage and Stability:	-20°C/1 year