

Antibodies Product Datasheet

Product Name:	EIF4E2 Antibody
Catalog Number:	C50569
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EIF4E2
Concentration:	1 mg/ml
Gene Name:	EIF4E2 EIF4EL3
Protein Name:	Eukaryotic translation initiation factor 4E type 2 (eIF-4E type 2) (eIF4E type 2) (Eukaryotic translation initiation factor 4E homologous protein) (Eukaryotic translation initiation factor 4E-like 3) (eIF4E-like protein 4E-LP) (mRNA cap-binding protein 4E
Other Name:	
Human Swiss Prot No.:	O60573
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EIF4E2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year