
Antibodies Product Datasheet

Product Name:	TFIIIC90 Antibody
Catalog Number:	C11871
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TF3C4. AA range:611-660
Concentration:	1 mg/ml
Gene Name:	GTF3C4
Protein Name:	General transcription factor 3C polypeptide 4
Other Name:	GTF3C4; General transcription factor 3C polypeptide 4; TF3C-delta; Transcription factor III C 90 kDa subunit; TFIIIC 90 kDa subunit; TFIIIC90; Transcription factor III C subunit delta
Human Swiss Prot No.:	Q9UKN8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TFIIIC90 Polyclonal Antibody detects endogenous levels of TFIIIC90 protein.
Storage and Stability:	-20°C/1 year