

Antibodies Product Datasheet

Product Name:	T2R39 Antibody
Catalog Number:	G756
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TAS2R39. AA range:218-267
Concentration:	1 mg/ml
Gene Name:	TAS2R39
Protein Name:	Taste receptor type 2 member 39
Other Name:	TAS2R39; Taste receptor type 2 member 39; T2R39; Taste receptor type 2 member 57; T2R57
Human Swiss Prot No.:	P59534
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	T2R39 Polyclonal Antibody detects endogenous levels of T2R39 protein.
Storage and Stability:	-20°C/1 year