

Antibodies Product Datasheet

Product Name:	Cleaved-Factor XIIIa (G39) Antibody
Catalog Number:	L0206
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FA13A. AA range:20-69
Concentration:	1 mg/ml
Gene Name:	F13A1
Protein Name:	Coagulation factor XIII A chain
Other Name:	F13A1; F13A; Coagulation factor XIII A chain; Coagulation factor XIIIa; Protein-glutamine gamma-glutamyltransferase A chain; Transglutaminase A chain
Human Swiss Prot No.:	P00488
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-Factor XIIIa (G39) Polyclonal Antibody detects endogenous levels of fragment of activated Factor XIIIa protein resulting from cleavage adjacent to G39.
Storage and Stability:	-20°C/1 year