

## Antibodies Product Datasheet

<b>Product Name:</b>	FTSJD2 Antibody
<b>Catalog Number:</b>	C50974
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal, Rabbit,IgG
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human FTSJD2
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	FTSJD2 KIAA0082 MTR1
<b>Protein Name:</b>	Cap-specific mRNA (nucleoside-2'-O-)-methyltransferase 1 (EC 2.1.1.57) (Cap1 2'O-ribose methyltransferase 1) (MTr1) (hMTr1) (FtsJ methyltransferase domain-containing protein 2) (Interferon-stimulated gene 95 kDa protein) (ISG95)
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	Q8N1G2
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of FTSJD2 at Human, Mouse,Rat
<b>Storage and Stability:</b>	Store at -20° C/1 year