

Antibodies Product Datasheet

Product Name:	Daxx Monoclonal Antibody
Catalog Number:	N0192
Reactivity:	Human
Applications:	WB,IF/ICC,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Daxx expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	DAXX
Protein Name:	Death domain-associated protein 6
Other Name:	DAXX; BING2; DAP6; Death domain-associated protein 6; Daxx; hDaxx; ETS1-associated protein 1; EAP1; Fas death domain-associated protein
Human Swiss Prot No.:	Q9UER7
Formulation:	Antibody are purified by protein G affinity chromatography. Liquid in PBS containing 50% glycerol and 0.03% sodium azide.
Specificity:	Daxx Monoclonal Antibody detects endogenous levels of Daxx protein.
Storage and Stability:	-20°C/1 year