

Antibodies Product Datasheet

Product Name:	TRIP13 Antibody
Catalog Number:	C11092
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TRIP13. AA range:383-432
Concentration:	1 mg/ml
Gene Name:	TRIP13
Protein Name:	Pachytene checkpoint protein 2 homolog
Other Name:	TRIP13; PCH2; Pachytene checkpoint protein 2 homolog; Human papillomavirus type 16 E1 protein-binding protein; 16E1-BP; HPV16 E1 protein-binding protein; Thyroid hormone receptor interactor 13; Thyroid receptor-interacting protein 13; TR-in
Human Swiss Prot No.:	Q15645
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TRIP13 Polyclonal Antibody detects endogenous levels of TRIP13 protein.
Storage and Stability:	-20°C/1 year