

Antibodies Product Datasheet

Product Name:	Cyclin E2 (Phospho-Thr392) Antibody
Catalog Number:	A0880
Reactivity:	Human, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cyclin E2 around the phosphorylation site of Thr392. AA range:355-404
Concentration:	1 mg/ml
Gene Name:	CCNE2
Protein Name:	G1/S-specific cyclin-E2
Other Name:	CCNE2; G1/S-specific cyclin-E2
Human Swiss Prot No.:	O96020
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Cyclin E2 (T392) Polyclonal Antibody detects endogenous levels of Cyclin E2 protein only when phosphorylated at T392.
Storage and Stability:	-20°C/1 year