

Antibodies Product Datasheet

Product Name:	DEFA5 Antibody
Catalog Number:	C49821
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DEFA5
Concentration:	1 mg/ml
Gene Name:	DEFA5 DEF5
Protein Name:	Defensin-5 (Defensin, alpha 5) (HD5(20-94)) [Cleaved into: HD5(23-94)HD5(29-94)HD5(56-94)HD5(63-94)]
Other Name:	
Human Swiss Prot No.:	Q01523
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DEFA5 at Human,Rat
Storage and Stability:	Store at -20° C/1 year