

Antibodies Product Datasheet

Product Name:	DYRK1A/B (Phospho-Tyr321/273) Antibody
Catalog Number:	A8093
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DYRK1A/B (Phospho-Tyr321/273)
Concentration:	1 mg/ml
Gene Name:	DYRK1A DYRK MNB MNBH
Protein Name:	DYRK1A/B (Phospho-Tyr321/273)
Other Name:	Dual specificity tyrosine-phosphorylation-regulated kinase 1A (EC 2.7.12.1) (Dual specificity YAK1-related kinase) (HP86) (Protein kinase minibrain homolog) (MNBH) (hMNB)
Human Swiss Prot No.:	Q13627
Formulation:	Liquid in PBS containing 50% glycerol, and 0.91% sodium azide.
Specificity:	This antibody detects endogenous levels of DYRK1A/B (Phospho-Tyr321/273) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year