

Antibodies Product Datasheet

Product Name:	KCNG3 Antibody
Catalog Number:	C17811
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KCNG3. AA range:183-232
Concentration:	1 mg/ml
Gene Name:	KCNG3
Protein Name:	Potassium voltage-gated channel subfamily G member 3
Other Name:	KCNG3; Potassium voltage-gated channel subfamily G member 3; Voltage-gated potassium channel subunit Kv10.1; Voltage-gated potassium channel subunit Kv6.3
Human Swiss Prot No.:	Q8TAE7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KCNG3 Polyclonal Antibody detects endogenous levels of KCNG3 protein.
Storage and Stability:	-20°C/1 year