

Antibodies Product Datasheet

Product Name:	CPB2 Antibody
Catalog Number:	R12-2647
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	CPB2 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CPB2.
Concentration:	1 mg/ml
Gene Name:	CPB2
Other Name:	Carboxypeptidase B2;Carboxypeptidase U;CPU;Plasma carboxypeptidase B;pCPB;Thrombin-activable fibrinolysis inhibitor;TAFI ;CPB2
Human Swiss Prot No.:	Q96IY4
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	CPB2 Antibodyantibody detects endogenous levels of CPB2.
Storage and Stability:	Stable at -20°C for at least 1 year.