

## Antibodies Product Datasheet

<b>Product Name:</b>	Anti-Human Calcitonin (ABT460) Mouse Monoclonal Antibody
<b>Catalog Number:</b>	V0204
<b>Clone ID:</b>	ABT460
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC
<b>Application Dilutions</b>	IHC 1:100-200
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>UniProt:</b>	CALCA
<b>Gene Name:</b>	P01258
<b>Other Name:</b>	Alpha CGRP, alpha type CGRP, CALC 1, CALC A, CALC1, calcitonin 1, calcitonin gene related peptide 1, calcitonin gene related peptide I, calcitonin related polypeptide alpha, CGRP 1, CGRP, CGRP I, CGRP1, CT, katalcalcin, KC
<b>Isotype</b>	IgG1, kappa
<b>Localization:</b>	Thyroid/ Medullarythyroidcarcinoma
<b>Positive Tissue Control:</b>	Cytoplasmic
<b>Immunogen</b>	Synthetic peptide
<b>Purification</b>	Immunogen affinity purified
<b>Formulation:</b>	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
<b>Specificity:</b>	The antibody can specifically recognize human Calcitonin protein.
<b>Storage and Stability:</b>	-20°C/1 year