
Antibodies Product Datasheet

Product Name:	AP-1 Monoclonal Antibody
Catalog Number:	N0028
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human AP-1 expressed in E. Coli.
Concentration:	
Gene Name:	JUN
Protein Name:	Transcription factor AP-1
Other Name:	JUN; Transcription factor AP-1; Activator protein 1; AP1; Proto-oncogene c-Jun; V-jun avian sarcoma virus 17 oncogene homolog; p39
Human Swiss Prot No.:	P05412
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	AP-1 Monoclonal Antibody detects endogenous levels of AP-1 protein.
Storage and Stability:	-20°C/1 year