

Antibodies Product Datasheet

Product Name:	ERp72 Antibody
Catalog Number:	C48073
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ERp72
Concentration:	1 mg/ml
Gene Name:	PDIA4 ERP70 ERP72
Protein Name:	Protein disulfide-isomerase A4 (EC 5.3.4.1) (Endoplasmic reticulum resident protein 70) (ER protein 70) (ERp70) (Endoplasmic reticulum resident protein 72) (ER protein 72) (ERp-72) (ERp72)
Other Name:	
Human Swiss Prot No.:	P13667
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ERp72 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year