

Antibodies Product Datasheet

Product Name:	Caspase 10 (p12, Cleaved-Ala416) Antibody
Catalog Number:	L0477
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Caspase 10 (p12, Cleaved-Ala416)
Concentration:	1 mg/ml
Gene Name:	CASP10 MCH4
Protein Name:	Caspase10
Other Name:	Caspase-10 (CASP-10) (EC 3.4.22.63) (Apoptotic protease Mch-4) (FAS-associated death domain protein interleukin-1B-converting enzyme 2) (FLICE2) (ICE-like apoptotic protease 4) [Cleaved into: Caspase-10 subunit p23/17; Caspase-10 subunit p12]
Human Swiss Prot No.:	Q92851
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Caspase 10 (p12, Cleaved-Ala416, protein was cleaved amino acid sequence between 415-416)
Storage and Stability:	-20°C/1 year