

Antibodies Product Datasheet

Product Name:	Anti-Human p120 Catenin (ABT361) Mouse Monoclonal Antibody
Catalog Number:	V0425
Clone ID:	ABT361
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	CTNND1
Gene Name:	O60716
Other Name:	Cadherin associated Src substrate, Cadherin-associated Src substrate, CAS, Catenin (cadherin associated protein) delta 1, Catenin delta 1, Catenin delta, Catenin delta-1, CTND1_HUMAN, CTNND 1, CTNND, CTNND1, delta 1 Catenin, KIAA0384, p120, P120 CAS, p120 catenin, P120 CTN, p120(cas), p120(ctn), P120CAS, P120CTN
Isotype	IgG2b, Kappa
Localization:	Breast/ Tonsil
Positive Tissue Control:	Cytoplasmic, Membranous
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human p120 catenin protein. In western blotting of Hela, A431 and HEK293 cell lysates, the antibody can label a band with molecular weight around 100~110 kDa.
Storage and Stability:	-20°C/1 year