

Antibodies Product Datasheet

Product Name:	TCP-1 θ Antibody
Catalog Number:	C19087
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CCT8. AA range:361-410
Concentration:	1 mg/ml
Gene Name:	CCT8
Protein Name:	T-complex protein 1 subunit theta
Other Name:	CCT8; C21orf112; CCTQ; KIAA0002; T-complex protein 1 subunit theta; TCP-1-theta; CCT-theta; Renal carcinoma antigen NY-REN-15
Human Swiss Prot No.:	P50990
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TCP-1 θ Polyclonal Antibody detects endogenous levels of TCP-1 θ protein.
Storage and Stability:	-20°C/1 year