

Antibodies Product Datasheet

Product Name:	IL-10R α Antibody
Catalog Number:	B1056
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IL-10R alpha. AA range:462-511
Concentration:	1 mg/ml
Gene Name:	IL10RA
Protein Name:	Interleukin-10 receptor subunit alpha
Other Name:	IL10RA; IL10R; Interleukin-10 receptor subunit alpha; IL-10 receptor subunit alpha; IL-10R subunit alpha; IL-10RA; CDw210a; Interleukin-10 receptor subunit 1; IL-10R subunit 1; IL-10R1; CD antigen CD210
Human Swiss Prot No.:	Q13651
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-10R α Polyclonal Antibody detects endogenous levels of IL-10R α protein.
Storage and Stability:	-20°C/1 year