

Antibodies Product Datasheet

Product Name:	Rab11-FIP3 Antibody
Catalog Number:	C18121
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RAB11FIP3. AA range:569-618
Concentration:	1 mg/ml
Gene Name:	RAB11FIP3
Protein Name:	Rab11 family-interacting protein 3
Other Name:	RAB11FIP3; ARFO1; KIAA0665; Rab11 family-interacting protein 3; FIP3-Rab11; Rab11-FIP3; Arfophilin-1; EF hands-containing Rab-interacting protein; Eferin; MU-MB-17.148
Human Swiss Prot No.:	O75154
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Rab11-FIP3 Polyclonal Antibody detects endogenous levels of Rab11-FIP3 protein.
Storage and Stability:	-20°C/1 year