
Antibodies Product Datasheet

Product Name:	Phospho-Stat3 (Y705) mouse Monoclonal Antibody(1E9)
Catalog Number:	N1623
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	WB 500-2000 IHC-p 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Phospho-Stat3 (Y705) at AA range of 640-720
Concentration:	1 mg/ml
Gene Name:	STAT3
Protein Name:	STAT3
Other Name:	STAT3
Human Swiss Prot No.:	P40763
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Stat3 (Y705) protein detects endogenous levels of STAT3
Storage and Stability:	-20°C/1 year