

Antibodies Product Datasheet

Product Name:	NPL Antibody
Catalog Number:	C50679
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human NPL
Concentration:	1 mg/ml
Gene Name:	NPL C1orf13
Protein Name:	N-acetylneuraminate lyase (NALase) (EC 4.1.3.3) (N-acetylneuraminate pyruvate-lyase) (N-acetylneuraminic acid aldolase) (Sialate lyase) (Sialate-pyruvate lyase) (Sialic acid aldolase) (Sialic acid lyase)
Other Name:	
Human Swiss Prot No.:	Q9BXD5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of NPL at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year