

Antibodies Product Datasheet

Product Name:	Anti-Human Cytokeratin 19&10 (ABT505) Mouse Monoclonal Antibody
Catalog Number:	V0304
Clone ID:	ABT505
Reactivity:	Human
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	KRT19
Gene Name:	P08727
Other Name:	40 kDa keratin intermediate filament, CK 19, CK-19, CK19, Cytokeratin 19, Cytokeratin-19, K19, K1C19_HUMAN, K1CS, Keratin 19, Keratin type I 40 kD, Keratin type I 40kD, Keratin type I cytoskeletal 19, Keratin, type I cytoskeletal 19, Keratin, type I, 40 kd, Keratin-19, KRT19, MGC15366
Isotype	IgG1, Kappa
Localization:	Kindeg/ Skin
Positive Tissue Control:	Cytoplasmic, Membranous
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can recognize human CK10 and CK19 protein, and shows no cross reaction with CK1, 5, 6, 7, 8, 13, 14, 15, 17, 18, 20.
Storage and Stability:	-20°C/1 year