

Antibodies Product Datasheet

Product Name:	Skp2 p45 Antibody
Catalog Number:	C0324
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SKP2/p45. AA range:366-415
Concentration:	1 mg/ml
Gene Name:	SKP2
Protein Name:	S-phase kinase-associated protein 2
Other Name:	SKP2; FBXL1; S-phase kinase-associated protein 2; Cyclin-A/CDK2-associated protein p45; F-box protein Skp2; F-box/LRR-repeat protein 1; p45skp2
Human Swiss Prot No.:	Q13309
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Skp2 p45 Polyclonal Antibody detects endogenous levels of Skp2 p45 protein.
Storage and Stability:	-20°C/1 year