

## Antibodies Product Datasheet

<b>Product Name:</b>	FGFR-4 Monoclonal Antibody
<b>Catalog Number:</b>	N0266
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB,IF/ICC,ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	Affinity purification
<b>Immunogen:</b>	Purified recombinant extracellular fragment of human FGFR-4 fused with hlgGfc tag expressed in HEK293 cell line.
<b>Concentration:</b>	
<b>Gene Name:</b>	FGFR4
<b>Protein Name:</b>	Fibroblast growth factor receptor 4
<b>Other Name:</b>	FGFR4; JTK2; TKF; Fibroblast growth factor receptor 4; FGFR-4; CD antigen CD334
<b>Human Swiss Prot No.:</b>	P22455
<b>Formulation:</b>	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
<b>Specificity:</b>	FGFR-4 Monoclonal Antibody detects endogenous levels of FGFR-4 protein.
<b>Storage and Stability:</b>	-20°C/1 year