

Antibodies Product Datasheet

Product Name:	UGT2B4 Antibody
Catalog Number:	C49391
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human UGT2B4
Concentration:	1 mg/ml
Gene Name:	UGT2B4 UGT2B11
Protein Name:	UDP-glucuronosyltransferase 2B4 (UDPGT 2B4) (EC 2.4.1.17) (HLUG25) (Hyodeoxycholic acid-specific UDPGT) (UDPGTh-1)
Other Name:	
Human Swiss Prot No.:	P06133
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of UGT2B4 at Human
Storage and Stability:	Store at -20° C/1 year