

## Antibodies Product Datasheet

<b>Product Name:</b>	Histone H4 (Phospho-Thr80) Antibody
<b>Catalog Number:</b>	A9730
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB: 1:1000-2000
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthetic Peptide of Histone H4 (Phospho Thr80)
<b>Concentration:</b>	
<b>Gene Name:</b>	HIST1H4A
<b>Protein Name:</b>	Histone H4
<b>Other Name:</b>	HIST1H4A; H4/A; H4FA; HIST1H4B; H4/I; H4FI; HIST1H4C; H4/G; H4FG; HIST1H4D; H4/B; H4FB; HIST1H4E; H4/J; H4FJ; HIST1H4F; H4/C; H4FC; HIST1H4H; H4/H; H4FH; HIST1H4I; H4/M; H4FM; HIST1H4J; H4/E; H4FE; HIST1H4K; H4/D; H4FD; HIST1H4L; H4/K; H4FK; HIST2H4A; H4/
<b>Human Swiss Prot No.:</b>	P62805
<b>Formulation:</b>	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Specificity:</b>	The antibody detects endogenous Histone H4 (Phospho Thr80) protein.
<b>Storage and Stability:</b>	-20°C/1 year