

Antibodies Product Datasheet

Product Name:	FRS2 Antibody
Catalog Number:	B8239
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FRS2. AA range:402-451
Concentration:	1 mg/ml
Gene Name:	FRS2
Protein Name:	Fibroblast growth factor receptor substrate 2
Other Name:	FRS2; Fibroblast growth factor receptor substrate 2; FGFR substrate 2; FGFR-signaling adaptor SNT; Suc1-associated neurotrophic factor target 1; SNT-1
Human Swiss Prot No.:	Q8WU20
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FRS2 Polyclonal Antibody detects endogenous levels of FRS2 protein.
Storage and Stability:	-20°C/1 year