

Antibodies Product Datasheet

Product Name:	DGK- ι Antibody
Catalog Number:	C10197
Reactivity:	Human, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DGKI. AA range:991-1040
Concentration:	1 mg/ml
Gene Name:	DGKI
Protein Name:	Diacylglycerol kinase ι
Other Name:	DGKI; Diacylglycerol kinase ι ; DAG kinase ι ; Diglyceride kinase ι ; DGK- ι
Human Swiss Prot No.:	O75912
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DGK- ι Polyclonal Antibody detects endogenous levels of DGK- ι protein.
Storage and Stability:	-20°C/1 year