

Antibodies Product Datasheet

Product Name:	PZR (Phospho-Tyr241) Antibody
Catalog Number:	A9609
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human PZR (Tyr241)
Concentration:	1 mg/ml
Gene Name:	MPZL1 PZR UNQ849/PRO1787
Protein Name:	PZR (Tyr241)
Other Name:	Myelin protein zero-like protein 1 (Protein zero-related)
Human Swiss Prot No.:	O95297
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat PZR (phospho-Tyr241)
Storage and Stability:	-20°C/1 year