
Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | DLK1 Monoclonal Antibody |
| Catalog Number: | N0199 |
| Reactivity: | Human |
| Applications: | WB,IF/ICC,FCM,ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications. |
| Modification: | Unmodified/ Total |
| Source: | Monoclonal Mouse |
| Purification: | Affinity purification |
| Immunogen: | Purified recombinant fragment of human DLK1 expressed in E. Coli. |
| Concentration: | |
| Gene Name: | DLK1 |
| Protein Name: | Protein delta homolog 1 |
| Other Name: | DLK1; DLK; Protein delta homolog 1; DLK-1; pG2 |
| Human Swiss Prot No.: | P80370 |
| Formulation: | Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol. |
| Specificity: | DLK1 Monoclonal Antibody detects endogenous levels of DLK1 protein. |
| Storage and Stability: | -20°C/1 year |