

Antibodies Product Datasheet

Product Name:	Anti-Human Cytokeratin 20 (ABT316) Mouse Monoclonal Antibody
Catalog Number:	V0306
Clone ID:	ABT316
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	KRT20
Gene Name:	P35900
Other Name:	CD20, CK 20, CK-20, CK20, Cytokeratin-20, Cytokeratin20, K1C20_HUMAN, K20, KA20, Keratin 20, keratin 20, type I, keratin 21, rat, homolog of, Keratin, Keratin type I cytoskeletal 20, Keratin-20, Keratin20, KRT 20, KRT 21, KRT20, KRT21, MGC35423, OTTHUMP00000164518, Protein IT, type I cytoskeletal 20
Isotype	IgG1, kappa
Localization:	Colon/ Appendix
Positive Tissue Control:	Cytoplasmic, Membranous
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human CK20 protein, and shows no cross reaction with CK1, 5, 6, 7, 8, 13, 14, 15, 17, 18.
Storage and Stability:	-20°C/1 year