

Antibodies Product Datasheet

Product Name:	PEA3 Antibody
Catalog Number:	C10614
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ETV4. AA range:281-330
Concentration:	1 mg/ml
Gene Name:	ETV4
Protein Name:	ETS translocation variant 4
Other Name:	ETV4; E1AF; PEA3; ETS translocation variant 4; Adenovirus E1A enhancer-binding protein; E1A-F; Polyomavirus enhancer activator 3 homolog; Protein PEA3
Human Swiss Prot No.:	P43268
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PEA3 Polyclonal Antibody detects endogenous levels of PEA3 protein.
Storage and Stability:	-20°C/1 year