

Antibodies Product Datasheet

Product Name:	S6A12 Antibody
Catalog Number:	C46127
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	SLC6A12
Protein Name:	Sodium- and chloride-dependent betaine transporter (BGT-1) (Na ⁽⁺⁾ /Cl ⁽⁻⁾ betaine/GABA transporter) (Solute carrier family 6 member 12)
Other Name:	
Human Swiss Prot No.:	P48065
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	S6A12 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year