
Antibodies Product Datasheet

Product Name:	Ob-R (Phospho-Tyr1141) Antibody
Catalog Number:	A8265
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human Ob-R (phospho Tyr1141)
Concentration:	1 mg/ml
Gene Name:	LEPR OBR
Protein Name:	Leptin receptor
Other Name:	LEPR; DB; OBR; Leptin receptor; LEP-R; HuB219; OB receptor; OB-R; CD antigen CD295
Human Swiss Prot No.:	P48357
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Ob-R (Y1141) Polyclonal Antibody detects endogenous levels of Ob-R protein only when phosphorylated at Y1141.
Storage and Stability:	-20°C/1 year