

Antibodies Product Datasheet

Product Name:	SPA9 Antibody
Catalog Number:	C46376
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	SERPINA9 GCET1 SERPINA11 UNQ692/PRO1337
Protein Name:	Serpin A9 (Centerin) (Germinal center B-cell-expressed transcript 1 protein)
Other Name:	
Human Swiss Prot No.:	Q86WD7
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	SPA9 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year