

Antibodies Product Datasheet

Product Name:	EPHB3 (Phospho-Tyr608) Antibody
Catalog Number:	A9518
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EPHB3 (Phospho Tyr608)
Concentration:	1 mg/ml
Gene Name:	EPHB3 ETK2 HEK2 TYRO6
Protein Name:	EPHB3 (Phospho Tyr608)
Other Name:	Ephrin type-B receptor 3 (EC 2.7.10.1) (EPH-like tyrosine kinase 2) (EPH-like kinase 2) (Embryonic kinase 2) (EK2) (hEK2) (Tyrosine-protein kinase TYRO6)
Human Swiss Prot No.:	P54753
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat EPHB3 (Phospho Tyr608)
Storage and Stability:	-20°C/1 year