

## Antibodies Product Datasheet

<b>Product Name:</b>	DCPS Antibody
<b>Catalog Number:</b>	C51018
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal, Rabbit,IgG
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human DCPS
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	DCPS DCS1 HINT5 HSPC015
<b>Protein Name:</b>	m7GpppX diphosphatase (EC 3.6.1.59) (DCS-1) (Hint-related 7meGMP-directed hydrolase) (Histidine triad protein member 5) (HINT-5) (Scavenger mRNA-decapping enzyme DcpS)
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	Q96C86
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of DCPS at Human, Mouse,Rat
<b>Storage and Stability:</b>	Store at -20° C/1 year