

Antibodies Product Datasheet

Product Name:	POTE15 Antibody
Catalog Number:	C17977
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human A26B1. AA range:481-530
Concentration:	1 mg/ml
Gene Name:	POTEB
Protein Name:	POTE ankyrin domain family member B
Other Name:	POTEB; A26B1; POTE15; POTE ankyrin domain family member B; ANKRD26-like family B member 1; Prostate; ovary, testis-expressed protein on chromosome 15; POTE-15
Human Swiss Prot No.:	Q6S5H4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	POTE15 Polyclonal Antibody detects endogenous levels of POTE15 protein.
Storage and Stability:	-20°C/1 year