
Antibodies Product Datasheet

Product Name:	RPS6KC1 Antibody
Catalog Number:	R12-3465
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	RPS6KC1 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RPS6KC1.
Concentration:	1 mg/ml
Gene Name:	RPS6KC1
Other Name:	Ribosomal protein S6 kinase delta-1; S6K-delta-1; 52 kDa ribosomal protein S6 kinase; Ribosomal S6 kinase-like protein with two PSK domains 118 kDa protein; SPHK1-binding protein; RPS6KC1; RPK118 Q96S38
Human Swiss Prot No.:	
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	RPS6KC1 Antibody antibody detects endogenous levels of RPS6KC1.
Storage and Stability:	Stable at -20°C for at least 1 year.