

Antibodies Product Datasheet

Product Name:	MBOA5 Antibody
Catalog Number:	C51177
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C-terminal MBOA5
Concentration:	1 mg/ml
Gene Name:	LPCAT3 MBOAT5 OACT5
Protein Name:	Lysophospholipid acyltransferase 5 (LPLAT 5) (EC 2.3.1.-) (1-acylglycerophosphocholine O-acyltransferase) (EC 2.3.1.23) (1-acylglycerophosphoserine O-acyltransferase) (EC 2.3.1.n6) (Lysophosphatidylcholine acyltransferase) (LPCAT) (Lyso-PC acyltransferase
Other Name:	Lysophospholipid acyltransferase 5 (LPLAT 5) (EC 2.3.1.-) (1-acylglycerophosphocholine O-acyltransferase) (EC 2.3.1.23) (1-acylglycerophosphoserine O-acyltransferase) (EC 2.3.1.n6) (Lysophosphatidylcholine acyltransferase) (LPCAT) (Lyso-PC acyltransferase
Human Swiss Prot No.:	Q6P1A2
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MBOA5 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year