

Antibodies Product Datasheet

Product Name:	Anti-Human Cytokeratin 15 (ABT477) Mouse Monoclonal Antibody
Catalog Number:	V0294
Clone ID:	ABT477
Reactivity:	Human
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	KRT15
Gene Name:	P19012
Other Name:	AI528832, CK 15, CK-15, CK15, Cytokeratin-15, Cytokeratin15, K15, K1C15_HUMAN, K1CO, Ka15, Keratin 15, Keratin 15 basic, Keratin 15 beta, Keratin, Keratin complex 1 acidic gene 15, Keratin type I cytoskeletal 15, Keratin-15, Keratin15, KRT 15, Krt1-15, KRT15, KRTB, KRTL15, Type I cytoskeletal 15, Type I keratin Ka15
Isotype	IgG2a, Kappa
Localization:	Tonsil
Positive Tissue Control:	Cytoplasmic, Membranous
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human CK15 protein. In immunohistochemistry on formalin-fixed, paraffin-embedded tissue sections, the antibody specifically labels the squamous epithelial cells and carcinomas, and shows no reaction with glandular e
Storage and Stability:	-20°C/1 year