

Antibodies Product Datasheet

Product Name:	CRSP70 Antibody
Catalog Number:	C10228
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MED26. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	MED26
Protein Name:	Mediator of RNA polymerase II transcription subunit 26
Other Name:	MED26; ARC70; CRSP7; Mediator of RNA polymerase II transcription subunit 26; Activator-recruited cofactor 70 kDa component; ARC70; Cofactor required for Sp1 transcriptional activation subunit 7; CRSP complex subunit 7; Mediator complex subu
Human Swiss Prot No.:	O95402
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CRSP70 Polyclonal Antibody detects endogenous levels of CRSP70 protein.
Storage and Stability:	-20°C/1 year