

Antibodies Product Datasheet

Product Name:	Rad52 (Phospho-Tyr104) Antibody
Catalog Number:	A1119
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RAD52 around the phosphorylation site of Tyr104. AA range:70-119
Concentration:	1 mg/ml
Gene Name:	RAD52
Protein Name:	DNA repair protein RAD52 homolog
Other Name:	RAD52; DNA repair protein RAD52 homolog
Human Swiss Prot No.:	P43351
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Rad52 (Y104) Polyclonal Antibody detects endogenous levels of Rad52 protein only when phosphorylated at Y104.
Storage and Stability:	-20°C/1 year