

Antibodies Product Datasheet

Product Name:	p53 (Acetyl-Lys370) Antibody
Catalog Number:	D12171
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized Acetyl-peptide derived from human p53 around the Acetylation site of Lys370. AA range:331-380
Concentration:	1 mg/ml
Gene Name:	TP53
Protein Name:	Cellular tumor antigen p53
Other Name:	TP53; P53; Cellular tumor antigen p53; Antigen NY-CO-13; Phosphoprotein p53; Tumor suppressor p53
Human Swiss Prot No.:	P04637
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Acetyl-p53 (K370) Polyclonal Antibody detects endogenous levels of p53 protein only when acetylated at K370.
Storage and Stability:	-20°C/1 year