
Antibodies Product Datasheet

Product Name:	Connexin-40 Antibody
Catalog Number:	C30205
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human GJA5. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	GJA5
Protein Name:	Gap junction alpha-5 protein
Other Name:	GJA5; Gap junction alpha-5 protein; Connexin-40; Cx40
Human Swiss Prot No.:	P36382
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Connexin-40 Polyclonal Antibody detects endogenous levels of Connexin-40 protein.
Storage and Stability:	-20°C/1 year