
Antibodies Product Datasheet

Product Name:	VCP Antibody
Catalog Number:	R12-3676
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	VCP Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human VCP.
Concentration:	1 mg/ml
Gene Name:	VCP
Other Name:	15S Mg(2+)-ATPase p97 subunit; TER ATPase; TERA; Transitional endoplasmic reticulum ATPase; Valosin containing protein; VCP
Human Swiss Prot No.:	P55072
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	VCP Antibody antibody detects endogenous levels of VCP.
Storage and Stability:	Stable at -20°C for at least 1 year.