

## Antibodies Product Datasheet

<b>Product Name:</b>	PRX I Antibody
<b>Catalog Number:</b>	C30292
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from the Internal region of human PRDX1. AA range:31-80
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	PRDX1
<b>Protein Name:</b>	Peroxiredoxin-1
<b>Other Name:</b>	PRDX1; PAGA; PAGB; TDPX2; Peroxiredoxin-1; Natural killer cell-enhancing factor A; NKEF-A; Proliferation-associated gene protein; PAG; Thioredoxin peroxidase 2; Thioredoxin-dependent peroxide reductase 2
<b>Human Swiss Prot No.:</b>	Q06830
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	PRX I Polyclonal Antibody detects endogenous levels of PRX I protein.
<b>Storage and Stability:</b>	-20°C/1 year